

Aesculap Neurosurgery

Yasargil Aneurysm Clip System

Long Clips



...the brand you trust in



YASARGIL Long Clip

- Extra long clip blades for treatment of giant aneurysms
- Box lock guidance mechanism prevents scissoring
- The pyramidal design, on the inner portion of the jaw blades, provides an atraumatic distribution of pressure across the width of the neck of the aneurysm

Phynox alloy acc. ISO 5832-7



Art. No.	Blade Length	Max. Opening	Closing Force	
FE861K	30 mm	19 mm	190 g / 1.86 N	
FE862K	35 mm	21 mm	200 g / 1.96 N	
FE863K	40 mm	23 mm	200 g / 1.96 N	

Applying Forceps for Long Clip

- With active, controllable locking mechanism; the applier can be used with or without the lock being activated, whatever the surgeon desires
- Marked with "Phynox Long Clip" for easy identification



Art. No.	FE560K
----------	--------

B | BRAUN
SHARING EXPERTISE

Aesculap AG & Co. KG

Am Aesculap-Platz
78532 Tuttlingen
Germany

Phone +49 7461 95-0
Fax +49 7461 95-2600

www.aesculap.de

All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogues and pricelists and to take legal actions.